

SUPER CPAP BATTERY PACKS

FREQUENTLY ASKED QUESTIONS

Is a prescription required and do insurance companies pay for the CPAP Battery Pack?

A prescription is not required. CPAP Battery Packs are not considered a medically necessary item and therefore will most likely not be covered by insurance.

Do CPAP Battery Packs work with all CPAP/BiLevel machine models?

CPAP Battery Packs are compatible with most makes and models of CPAP and BiLevel machines; however, use of a sine wave inverter or DC converter may be required.

When is a sine wave inverter or DC converter required?

The use of a sine wave inverter or DC converter is required for all ResMed model machines or if the device does not have a DC input port.

How long will the Battery Pack run a CPAP machine?

Run times vary depending upon several factors including device type, pressure settings and altitude. Please note the C-100 Battery Pack can be used in a single battery configuration or a pigtailed double battery configuration.

Listed run times per charge based on Respironics System One and ResMed S-9 models with no humidification:

C-100 (single configuration): 20-25 hours C-100 (double configuration): 40-45 hours

C-150: 15-20 hours C-222: 25-30 hours C-444: 65+ hours

Use of a sine wave inverter, DC converter or a heated humidifier will reduce the expected run time for all battery models.

How long does it take to recharge the Battery Pack?

CPAP Battery Packs are made from high capacity lithium-ion cells so there is not the typical “memory loss” issues that are associated with other types of batteries. The battery does not need to be fully discharged before it is recharged. Battery recharge times from full discharge to full recharge:

C-100 (per battery): 3-4 hours C-150: 4-5 hours C-222: 5-6 hours C-444: 6-7 hours

Can the Battery Pack be recharged via the cigarette lighter socket of a vehicle?

Yes. Battery Packs can be recharged via the cigarette lighter socket of a vehicle by using a sine wave inverter or DC converter in conjunction with the Battery Pack.

Can CPAP Battery Packs be used to power other devices?

YES. Super CPAP Battery Packs can also be used to power other devices that operate from 12V DC current such as portable DVD players, iPods, laptops, iPads, cell phones, and more! Use of a sine wave inverter or DC converter may be required.

Where can I purchase Super CPAP Battery Packs?

Battery Power Solutions is a business-to-business company and does not sell directly to consumers. Please visit our website at www.batterypowersolutions.net and click on the “Where to Buy” tab to locate an authorized Super CPAP Battery Pack dealer near you.

1-877-445-5228
info@batterypowersolutions.net

